

Antibodies Product Datasheet

Product Name:	SLC39A3 Antibody
Catalog Number:	C51041
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC39A3
Concentration:	1 mg/ml
Gene Name:	SLC39A3 ZIP3
Protein Name:	Zinc transporter ZIP3 (Solute carrier family 39 member 3) (Zrt- and Irt-like protein 3) (ZIP-3)
Other Name:	
Human Swiss Prot No.:	Q9BRY0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC39A3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year