

Antibodies Product Datasheet

Product Name:	MASP-2 Antibody
Catalog Number:	C16548
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MASP2. AA range:227-276
Concentration:	1 mg/ml
Gene Name:	MASP2
Protein Name:	Mannan-binding lectin serine protease 2
Other Name:	MASP2; Mannan-binding lectin serine protease 2; MBL-associated serine protease 2; Mannose-binding protein-associated serine protease 2; MASP-2
Human Swiss Prot No.:	O00187
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MASP-2 Polyclonal Antibody detects endogenous levels of MASP-2 protein.
Storage and Stability:	-20°C/1 year