
Antibodies Product Datasheet

Product Name:	Blk Monoclonal Antibody
Catalog Number:	N0063
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Blk expressed in E. Coli.
Concentration:	
Gene Name:	BLK
Protein Name:	Tyrosine-protein kinase Blk
Other Name:	BLK; Tyrosine-protein kinase Blk; B lymphocyte kinase; p55-Blk
Human Swiss Prot No.:	P51451
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Blk Monoclonal Antibody detects endogenous levels of Blk protein.
Storage and Stability:	-20°C/1 year