

Antibodies Product Datasheet

Product Name:	PECR Antibody
Catalog Number:	C45487
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	PECR PRO1004
Protein Name:	Peroxisomal trans-2-enoyl-CoA reductase (TERP) (EC 1.3.1.38) (2,4-dienoyl-CoA reductase-related protein) (DCR-RP) (HPDHase) (pVI-ARL)
Other Name:	
Human Swiss Prot No.:	Q9BY49
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	PECR Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year