

Antibodies Product Datasheet

Product Name:	c-Myb (Phospho-Ser532) Antibody
Catalog Number:	A8378
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Myb around the phosphorylation site of Ser532. AA range:496-545
Concentration:	1 mg/ml
Gene Name:	MYB
Protein Name:	Transcriptional activator Myb
Other Name:	MYB; Transcriptional activator Myb; Proto-oncogene c-Myb
Human Swiss Prot No.:	P10242
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-c-Myb (S532) Polyclonal Antibody detects endogenous levels of c-Myb protein only when phosphorylated at S532.
Storage and Stability:	-20°C/1 year