

Antibodies Product Datasheet

Product Name:	PR Monoclonal Antibody(Z15)
Catalog Number:	N1506
Reactivity:	Human,Mouse,Rat
Applications:	IHC-p,IF/ICC
Dilutions:	IHC: 1:200 IF: 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of PR
Concentration:	
Gene Name:	PGR
Protein Name:	Progesterone receptor
Other Name:	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3
Human Swiss Prot No.:	P06401
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous PR protein.
Storage and Stability:	-20°C/1 year