
Antibodies Product Datasheet

Product Name:	M3K13 Antibody
Catalog Number:	R12-2982
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	M3K13 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human M3K13.
Concentration:	1 mg/ml
Gene Name:	MAP3K13
Other Name:	Mitogen-activated protein kinase kinase kinase 13; Mixed lineage kinase; MLK; Leucine zipper-bearing kinase; MAP3K13; LZK
Human Swiss Prot No.:	O43283
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	M3K13 Antibody antibody detects endogenous levels of M3K13.
Storage and Stability:	Stable at -20°C for at least 1 year.