

Antibodies Product Datasheet

Product Name:	ValRS Antibody
Catalog Number:	C50982
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ValRS
Concentration:	1 mg/ml
Gene Name:	VAR5 G7A VAR52
Protein Name:	Valine--tRNA ligase (EC 6.1.1.9) (Protein G7a) (Valyl-tRNA synthetase) (ValRS)
Other Name:	
Human Swiss Prot No.:	P26640
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ValRS at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year