
Antibodies Product Datasheet

Product Name:	α skeletal muscle actin Monoclonal Antibody(4B11)
Catalog Number:	N1414
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICCIP
Dilutions:	WB: 1:500-10000 IP:1:200 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of α skeletal muscle actin
Concentration:	
Gene Name:	ACTA1
Protein Name:	Actin alpha skeletal muscle
Other Name:	ACTA1; ACTA; Actin, alpha skeletal muscle; Alpha-actin-1
Human Swiss Prot No.:	P68133
Formulation:	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous α Skeletal Muscle Actin protein.
Storage and Stability:	-20°C/1 year