

## Antibodies Product Datasheet

<b>Product Name:</b>	eIF3L Antibody
<b>Catalog Number:</b>	C11925
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human IF3EI. AA range:1-50
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	EIF3L
<b>Protein Name:</b>	Eukaryotic translation initiation factor 3 subunit L
<b>Other Name:</b>	EIF3L; EIF3EIP; EIF3S6IP; HSPC021; HSPC025; MSTP005; Eukaryotic translation initiation factor 3 subunit L; eIF3L; Eukaryotic translation initiation factor 3 subunit 6-interacting protein; Eukaryotic translation initiation factor 3 subunit E
<b>Human Swiss Prot No.:</b>	Q9Y262
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	eIF3L Polyclonal Antibody detects endogenous levels of eIF3L protein.
<b>Storage and Stability:</b>	-20°C/1 year