
Antibodies Product Datasheet

Product Name:	AMPK α 1 mouse Monoclonal Antibody(9G3)
Catalog Number:	N1660
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of AMPK α 1 of PRKAA1
Concentration:	1 mg/ml
Gene Name:	PRKAA1
Protein Name:	PRKAA1
Other Name:	PRKAA1
Human Swiss Prot No.:	Q13131
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AMPK α 1 protein detects endogenous levels of PRKAA1
Storage and Stability:	-20°C/1 year