

Antibodies Product Datasheet

Product Name:	SR-7 Antibody
Catalog Number:	G018
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HTR7. AA range:391-440
Concentration:	1 mg/ml
Gene Name:	HTR7
Protein Name:	5-hydroxytryptamine receptor 7
Other Name:	HTR7; 5-hydroxytryptamine receptor 7; 5-HT-7; 5-HT7; 5-HT-X; Serotonin receptor 7
Human Swiss Prot No.:	P34969
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SR-7 Polyclonal Antibody detects endogenous levels of SR-7 protein.
Storage and Stability:	-20°C/1 year