

Antibodies Product Datasheet

Product Name:	Beta-catenin(PT2250) mouse mAb
Catalog Number:	N1868
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG2b, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Beta-catenin
Concentration:	
Gene Name:	CTNNB1 CTNNB OK/SW-cl.35 PRO2286
Protein Name:	Catenin beta-1 (Beta-catenin)
Other Name:	
Human Swiss Prot No.:	P35222
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Catenin beta-1 (Beta-catenin). Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH9.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year