

Antibodies Product Datasheet

Product Name:	P53 (Mono Methyl Lys370) Monoclonal Antibody(Q59)
Catalog Number:	N1530
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p: 1:50-1:200 WB 500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of P53 (Mono Methyl Lys370)
Concentration:	1 mg/ml
Gene Name:	TP53
Protein Name:	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
Other Name:	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
Human Swiss Prot No.:	P04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous P53 (Mono Methyl Lys370) protein
Storage and Stability:	-20°C/1 year