

Antibodies Product Datasheet

Product Name:	Anti-Human CD1a (ABT291) Mouse Monoclonal Antibody
Catalog Number:	V0223
Clone ID:	ABT291
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CD1A
Gene Name:	P06126
Other Name:	CD 1a, CD1, CD1a, CD1A Antigen, CD1A antigen, a polypeptide, CD1a molecule, CD1A_HUMAN, cluster of differentiation 1 A, cortical thymocyte antigen CD1A, differentiation antigen CD1 alpha 3, epidermal dendritic cell marker CD1a, FCB 6, FCB6, HTA 1, HTA1, hTa1 thymocyte antigen, OTTHUMP00000018907, R 4, R4, T 6, T-cell surface antigen T6/Leu-6, T-cell surface glycoprotein CD1a, T6
Isotype	IgG2b, Kappa
Localization:	Tonsil/ skin
Positive Tissue Control:	Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CD1a protein.
Storage and Stability:	-20°C/1 year