

Antibodies Product Datasheet

Product Name:	HMGCL Antibody
Catalog Number:	C51167
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal HMGCL
Concentration:	1 mg/ml
Gene Name:	HMGCL
Protein Name:	Hydroxymethylglutaryl-CoA lyase, mitochondrial (HL) (HMG-CoA lyase) (EC 4.1.3.4) (3-hydroxy-3-methylglutarate-CoA lyase)
Other Name:	Hydroxymethylglutaryl-CoA lyase, mitochondrial (HL) (HMG-CoA lyase) (EC 4.1.3.4) (3-hydroxy-3-methylglutarate-CoA lyase)
Human Swiss Prot No.:	P35914
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HMGCL at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year