

Antibodies Product Datasheet

Product Name:	TIS11D Antibody
Catalog Number:	C10632
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TISD. AA range:207-256
Concentration:	1 mg/ml
Gene Name:	ZFP36L2
Protein Name:	Zinc finger protein 36 C3H1 type-like 2
Other Name:	ZFP36L2; BRF2; ERF2; RNF162C; TIS11D; Zinc finger protein 36; C3H1 type-like 2; ZFP36-like 2; Butyrate response factor 2; EGF-response factor 2; ERF-2; Protein TIS11D
Human Swiss Prot No.:	P47974
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TIS11D Polyclonal Antibody detects endogenous levels of TIS11D protein.
Storage and Stability:	-20°C/1 year