

Antibodies Product Datasheet

Product Name:	UGCGL1 Antibody
Catalog Number:	C48473
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human UGCGL1
Concentration:	1 mg/ml
Gene Name:	UGGT1 GT UGCGL1 UGGT UGT1 UGTR
Protein Name:	UDP-glucose:glycoprotein glucosyltransferase 1 (UGT1) (hUGT1) (EC 2.4.1.-) (UDP--Glc:glycoprotein glucosyltransferase) (UDP-glucose ceramide glucosyltransferase-like 1)
Other Name:	
Human Swiss Prot No.:	Q9NYU2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of UGCGL1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year