

## Antibodies Product Datasheet

<b>Product Name:</b>	Flg Antibody
<b>Catalog Number:</b>	B0481
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human FGFR1. AA range:626-675
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	FGFR1
<b>Protein Name:</b>	Fibroblast growth factor receptor 1
<b>Other Name:</b>	FGFR1; BFGFR; CEK; FGFBR; FLG; FLT2; HBGFR; Fibroblast growth factor receptor 1; FGFR-1; Basic fibroblast growth factor receptor 1; BFGFR; bFGF-R-1; Fms-like tyrosine kinase 2; FLT-2; N-sam; Proto-oncogene c-Fgr; CD antigen CD331
<b>Human Swiss Prot No.:</b>	P11362
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Flg Polyclonal Antibody detects endogenous levels of Flg protein.
<b>Storage and Stability:</b>	-20°C/1 year