

Antibodies Product Datasheet

Product Name:	BBOX1 Antibody
Catalog Number:	C50344
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BBOX1
Concentration:	1 mg/ml
Gene Name:	BBOX1 BBH BBOX
Protein Name:	Gamma-butyrobetaine dioxygenase (EC 1.14.11.1) (Gamma-butyrobetaine hydroxylase) (Gamma-BBH) (Gamma-butyrobetaine,2-oxoglutarate dioxygenase)
Other Name:	
Human Swiss Prot No.:	O75936
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BBOX1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year