

Antibodies Product Datasheet

Product Name:	MEK2 Monoclonal Antibody(2C3)
Catalog Number:	N1537
Reactivity:	Human,Rat,Mouse
Applications:	WB
Dilutions:	Western Blot: 1/500 - 1/2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MEK-2 around the non-phosphorylation site of T394.
Concentration:	1 mg/ml
Gene Name:	MAP2K2
Protein Name:	Dual specificity mitogen-activated protein kinase kinase 2 (MAP kinase kinase 2) (MAPKK 2) (EC 2.7.12.2) (ERK activator kinase 2) (MAPK/ERK kinase 2) (MEK 2)
Other Name:	MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity mitogen-activated protein kinase kinase 2; MAP kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2; MEK 2
Human Swiss Prot No.:	P36507
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MEK2 protein detects endogenous levels of MEK2
Storage and Stability:	-20°C/1 year