

Antibodies Product Datasheet

Product Name:	Acetyl P53(K382) mouse Monoclonal Antibody(5H10)
Catalog Number:	N1702
Reactivity:	Human,Mouse,Rat
Applications:	IHC-p,IF(paraffin section)
Dilutions:	WB 500-2000 IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of human Acetyl P53(K382)
Concentration:	1 mg/ml
Gene Name:	TP53 P53
Protein Name:	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
Other Name:	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
Human Swiss Prot No.:	P04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Acetyl P53(K382) protein in human, Acetyl P53(K379) protein in mouse, Acetyl P53(K380) protein in Rat
Storage and Stability:	-20°C/1 year