
Antibodies Product Datasheet

Product Name:	FLT3 Antibody
Catalog Number:	R12-2770
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	FLT3 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FLT3.
Concentration:	1 mg/ml
Gene Name:	FLT3
Other Name:	Tyrosine-protein kinase receptor FLT3; Stem cell tyrosine kinase 1; STK-1; CD135 antigen
Human Swiss Prot No.:	P36888
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	FLT3 Antibody antibody detects endogenous levels of FLT3.
Storage and Stability:	Stable at -20°C for at least 1 year.