

Antibodies Product Datasheet

Product Name:	PLK1 Mouse mAb(Mix-mA)
Catalog Number:	N1495
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PLK1
Concentration:	1 mg/ml
Gene Name:	PLK1 PLK
Protein Name:	PLK1
Other Name:	Serine/threonine-protein kinase PLK1 (EC 2.7.11.21) (Polo-like kinase 1) (PLK-1) (Serine/threonine-protein kinase 13) (STPK13)
Human Swiss Prot No.:	P53350
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.68% sodium azide.
Specificity:	This antibody detects endogenous levels of PLK1 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year