

Antibodies Product Datasheet

Product Name:	FOS Antibody
Catalog Number:	R12-2772
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	FOS Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOS.
Concentration:	1 mg/ml
Gene Name:	FOS
Other Name:	Cellular oncogene fos; FOS; G0/G1 switch regulatory protein 7; G0S7; Proto-oncogene protein c-fos
Human Swiss Prot No.:	P01100
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	FOS Antibodyantibody detects endogenous levels of FOS.
Storage and Stability:	Stable at -20°C for at least 1 year.