

Antibodies Product Datasheet

Product Name:	SPAM1 Antibody
Catalog Number:	C50307
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SPAM1
Concentration:	1 mg/ml
Gene Name:	SPAM1 HYAL3 PH20
Protein Name:	Hyaluronidase PH-20 (Hyal-PH20) (EC 3.2.1.35) (Hyaluronoglucosaminidase PH-20) (Sperm adhesion molecule 1) (Sperm surface protein PH-20)
Other Name:	
Human Swiss Prot No.:	P38567
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SPAM1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year