

Antibodies Product Datasheet

Product Name:	c-Kit Monoclonal Antibody
Catalog Number:	N0155
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of c-Kit expressed in E. Coli.
Concentration:	
Gene Name:	KIT
Protein Name:	Mast/stem cell growth factor receptor
Other Name:	KIT; SCFR; Mast/stem cell growth factor receptor Kit; SCFR; Piebald trait protein; PBT; Proto-oncogene c-Kit; Tyrosine-protein kinase Kit; p145 c-kit; v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog; CD antigen CD117
Human Swiss Prot No.:	P10721
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	c-Kit Monoclonal Antibody detects endogenous levels of c-Kit protein.
Storage and Stability:	-20°C/1 year