

Antibodies Product Datasheet

Product Name:	PS-1 Antibody
Catalog Number:	C0308
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Presenilin 1. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	PSEN1
Protein Name:	Presenilin-1
Other Name:	PSEN1; AD3; PS1; PSNL1; Presenilin-1; PS-1; Protein S182
Human Swiss Prot No.:	P49768
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PS-1 Polyclonal Antibody detects endogenous levels of PS-1 protein.
Storage and Stability:	-20°C/1 year