

Antibodies Product Datasheet

Product Name:	FDPS Antibody
Catalog Number:	C48149
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FDPS
Concentration:	1 mg/ml
Gene Name:	FDPS FPS KIAA1293
Protein Name:	Farnesyl pyrophosphate synthase (FPP synthase) (FPS) (EC 2.5.1.10) ((2E,6E)-farnesyl diphosphate synthase) (Dimethylallyltranstransferase) (EC 2.5.1.1) (Farnesyl diphosphate synthase) (Geranyltranstransferase)
Other Name:	
Human Swiss Prot No.:	P14324
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FDPS at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year