
Antibodies Product Datasheet

Product Name:	Raf-B Monoclonal Antibody
Catalog Number:	N0547
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Raf-B expressed in E. Coli.
Concentration:	
Gene Name:	BRAF
Protein Name:	Serine/threonine-protein kinase B-raf
Other Name:	BRAF; BRAF1; RAFB1; Serine/threonine-protein kinase B-raf; Proto-oncogene B-Raf; p94; v-Raf murine sarcoma viral oncogene homolog B1
Human Swiss Prot No.:	P15056
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Raf-B Monoclonal Antibody detects endogenous levels of Raf-B protein.
Storage and Stability:	-20°C/1 year