

Antibodies Product Datasheet

Product Name:	NFκB-p52 (Cleaved-A454) Antibody
Catalog Number:	L0474
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NFκB-p52 (Cleaved-A454)
Concentration:	1 mg/ml
Gene Name:	NFKB2 LYT10
Protein Name:	NFκB-p52 (Cleaved-A454)
Other Name:	Nuclear factor NF-kappa-B p100 subunit (DNA-binding factor KBF2) (H2TF1) (Lymphocyte translocation chromosome 10 protein) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2) (Oncogene Lyt-10) (Lyt10) [Cleaved into: Nuclear factor NF-kap
Human Swiss Prot No.:	Q00653
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Cleaved-NFκB-p52 (A454, protein was cleaved amino acid sequence between 454-455)
Storage and Stability:	-20°C/1 year