

## Antibodies Product Datasheet

<b>Product Name:</b>	STOML1 Antibody
<b>Catalog Number:</b>	C50900
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human STOML1
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	STOML1 SLP1 UNC24 MSTP019
<b>Protein Name:</b>	Stomatin-like protein 1 (SLP-1) (EPB72-like protein 1) (Protein unc-24 homolog) (Stomatin-related protein) (STORP)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q9UBI4
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of STOML1 at Human, Mouse
<b>Storage and Stability:</b>	Store at -20° C/1 year