

Antibodies Product Datasheet

Product Name:	PEG10 Monoclonal Antibody
Catalog Number:	N0514
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of PEG10 (aa1-120) expressed in E. Coli.
Concentration:	
Gene Name:	PEG10
Protein Name:	Retrotransposon-derived protein PEG10
Other Name:	PEG10; EDR; KIAA1051; MAR2; MART2; MEF3L1; RGAG3; Retrotransposon-derived protein PEG10; Embryonal carcinoma differentiation-regulated protein; Mammalian retrotransposon-derived protein 2; Myelin expression factor 3-like protein 1; MEF3-lik
Human Swiss Prot No.:	Q86TG7
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PEG10 Monoclonal Antibody detects endogenous levels of PEG10 protein.
Storage and Stability:	-20°C/1 year