

Antibodies Product Datasheet

Product Name:	KAT6B Antibody
Catalog Number:	C44349
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1170-1250
Concentration:	1 mg/ml
Gene Name:	KAT6B KIAA0383 MORF MOZ2 MYST4
Protein Name:	Histone acetyltransferase KAT6B (EC 2.3.1.48) (Histone acetyltransferase MOZ2) (MOZ, YBF2/SAS3, SAS2 and TIP60 protein 4) (MYST-4) (Monocytic leukemia zinc finger protein-related factor)
Other Name:	
Human Swiss Prot No.:	Q8WYB5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KAT6B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year