

Antibodies Product Datasheet

Product Name:	VPAC1 Antibody
Catalog Number:	C30281
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human VIPR1. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	VIPR1
Protein Name:	Vasoactive intestinal polypeptide receptor 1
Other Name:	VIPR1; Vasoactive intestinal polypeptide receptor 1; VIP-R-1; Pituitary adenylate cyclase-activating polypeptide type II receptor; PACAP type II receptor; PACAP-R-2; PACAP-R2; VPAC1
Human Swiss Prot No.:	P32241
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	VPAC1 Polyclonal Antibody detects endogenous levels of VPAC1 protein.
Storage and Stability:	-20°C/1 year