

Antibodies Product Datasheet

Product Name:	Caspase-1 Antibody
Catalog Number:	C30122
Reactivity:	Human, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 WB 1:500-2000, IHC 1:50-300, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human CASP1. AA range:350-400
Concentration:	1 mg/ml
Gene Name:	CASP1 IL1BC IL1BCE
Protein Name:	Caspase1
Other Name:	caspase 1, apoptosis-related cysteine peptidase (interleukin 1, beta, convertase)
Human Swiss Prot No.:	P29466
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Caspase-1 Polyclonal Antibody detects endogenous levels of Caspase-1
Storage and Stability:	-20°C/1 year