

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Rabbit Anti Chicken IgY (H&L) Biotin
<b>Catalog Number:</b>	SA0363
<b>Reactivity:</b>	Chicken
<b>Conjugate:</b>	Biotin
<b>Applications:</b>	ELISA, WB, IHC, ICC
<b>Dilutions:</b>	ELISA (1/5000 - 1/20000), WB (1/5000 - 1/20000), IHC (1/500 - 1/1000), ICC (1/500 - 1/1000)
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Concentration:</b>	1mg/mL
<b>Other Name:</b>	Rabbit Anti-Chicken IgY (H&L) Biotin
<b>Formulation:</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
<b>Storage and Stability:</b>	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing