

## Antibodies Product Datasheet

<b>Product Name:</b>	CYP24A1 Antibody
<b>Catalog Number:</b>	C12253
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 24A1. AA range:411-460
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CYP24A1
<b>Protein Name:</b>	1,25-dihydroxyvitamin D(3) 24-hydroxylase mitochondrial
<b>Other Name:</b>	CYP24A1; CYP24; 1; 25-dihydroxyvitamin D(3) 24-hydroxylase, mitochondrial; 24-OHase; Vitamin D(3) 24-hydroxylase; Cytochrome P450 24A1; Cytochrome P450-CC24
<b>Human Swiss Prot No.:</b>	Q07973
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	CYP24A1 Polyclonal Antibody detects endogenous levels of CYP24A1 protein.
<b>Storage and Stability:</b>	-20°C/1 year