

Antibodies Product Datasheet

Product Name:	CYP39A1 Antibody
Catalog Number:	C21061
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CYP39A1.
Concentration:	1 mg/ml
Gene Name:	CYP39A1
Protein Name:	24-hydroxycholesterol 7-alpha-hydroxylase
Other Name:	CYP39A1; 24-hydroxycholesterol 7-alpha-hydroxylase; Cytochrome P450 39A1; hCYP39A1; Oxysterol 7-alpha-hydroxylase
Human Swiss Prot No.:	Q9NYL5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP39A1 Polyclonal Antibody detects endogenous levels of CYP39A1 protein.
Storage and Stability:	-20°C/1 year