

Antibodies Product Datasheet

Product Name:	CKR-2 Antibody
Catalog Number:	G048
Reactivity:	Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from CKR-2 . at AA range: 250-330
Concentration:	1 mg/ml
Gene Name:	CCR2
Protein Name:	C-C chemokine receptor type 2
Other Name:	CCR2; CMKBR2; C-C chemokine receptor type 2; C-C CKR-2; CC-CKR-2; CCR-2; CCR2; Monocyte chemoattractant protein 1 receptor; MCP-1-R; CD antigen CD192
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CKR-2 Polyclonal Antibody detects endogenous levels of CKR-2 protein.
Storage and Stability:	-20°C/1 year