

Antibodies Product Datasheet

Product Name:	CHSTE Antibody
Catalog Number:	C43082
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CHST14 D4ST1 UNQ1925/PRO4400
Protein Name:	Carbohydrate sulfotransferase 14 (EC 2.8.2.35) (Dermatan 4-sulfotransferase 1) (D4ST-1) (hD4ST1)
Other Name:	
Human Swiss Prot No.:	Q8NCH0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CHSTE Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year