

Antibodies Product Datasheet

Product Name:	VE-Cadherin mouse Monoclonal Antibody(8E5)
Catalog Number:	N1655
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of VE-Cadherin at AA range of 670-750
Concentration:	1 mg/ml
Gene Name:	CDH5
Protein Name:	CDH5
Other Name:	CDH5
Human Swiss Prot No.:	P33151
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	VE-Cadherin protein detects endogenous levels of CDH5
Storage and Stability:	-20°C/1 year