

Antibodies Product Datasheet

Product Name:	TFIIH p44 Antibody
Catalog Number:	C11010
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TF2H2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	GTF2H2
Protein Name:	General transcription factor IIH subunit 2
Other Name:	GTF2H2; BTF2P44; General transcription factor IIH subunit 2; Basic transcription factor 2 44 kDa subunit; BTF2 p44; General transcription factor IIH polypeptide 2; TFIIH basal transcription factor complex p44 subunit
Human Swiss Prot No.:	Q13888
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TFIIH p44 Polyclonal Antibody detects endogenous levels of TFIIH p44 protein.
Storage and Stability:	-20°C/1 year