

Antibodies Product Datasheet

Product Name:	HptII mouse mAb
Catalog Number:	N0990
Reactivity:	Transgenic
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	
Concentration:	1 mg/ml
Gene Name:	hptii
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transgenic Only.
Storage and Stability:	-20°C/1 year