

## Antibodies Product Datasheet

<b>Product Name:</b>	$\alpha$ -tubulin Antibody
<b>Catalog Number:</b>	C40023
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB: 1:2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Recombinant Protein of Tubulin alpha-1A chain
<b>Concentration:</b>	
<b>Gene Name:</b>	TUBA1A
<b>Protein Name:</b>	Tubulin alpha-1A chain
<b>Other Name:</b>	YM3217
<b>Human Swiss Prot No.:</b>	Q71U36/P68363
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous $\alpha$ -tubulin protein.
<b>Storage and Stability:</b>	-20°C/1 year