

Antibodies Product Datasheet

Product Name:	IL-6R α Antibody
Catalog Number:	C30041
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL6R. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	IL6R
Protein Name:	Interleukin-6 receptor subunit alpha
Other Name:	IL6R; Interleukin-6 receptor subunit alpha; IL-6 receptor subunit alpha; IL-6R subunit alpha; IL-6R-alpha; IL-6RA; IL-6R 1; Membrane glycoprotein 80; gp80; CD126
Human Swiss Prot No.:	P08887
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-6R α Polyclonal Antibody detects endogenous levels of IL-6R α protein.
Storage and Stability:	-20°C/1 year