
Antibodies Product Datasheet

Product Name:	BIOTIN*Polyclonal Rabbit Anti Human IgG(h+l)
Catalog Number:	SA0297
Reactivity:	HumanIgG(h+l)
Conjugate:	Biotin
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	BIOTIN*Polyclonal Rabbit Anti-Human IgG(h+l)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year