
Antibodies Product Datasheet

Product Name:	Rad50 (Phospho-Ser635) Antibody
Catalog Number:	A8409
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Rad50 (Ser635)
Concentration:	1 mg/ml
Gene Name:	RAD50
Protein Name:	Rad50 (Ser635)
Other Name:	DNA repair protein RAD50 (hRAD50) (EC 3.6.-.-)
Human Swiss Prot No.:	Q92878
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Rad50 (phospho-Ser635)
Storage and Stability:	-20°C/1 year