

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Cytokeratin 18 Monoclonal Antibody(1B3)
<b>Catalog Number:</b>	N1360
<b>Reactivity:</b>	Human,Mouse,Rat
<b>Applications:</b>	WB,IF/ICC
<b>Dilutions:</b>	WB: 1:1000 -3000 IF: 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Recombinant Protein of Keratin type I cytoskeletal 18
<b>Concentration:</b>	
<b>Gene Name:</b>	KRT18
<b>Protein Name:</b>	Keratin type I cytoskeletal 18
<b>Other Name:</b>	KRT18; CYK18; PIG46; Keratin, type I cytoskeletal 18; Cell proliferation-inducing gene 46 protein; Cytokeratin-18; CK-18; Keratin-18; K18
<b>Human Swiss Prot No.:</b>	P05783
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous Cytokeratin 18 protein.
<b>Storage and Stability:</b>	-20°C/1 year