

Antibodies Product Datasheet

Product Name:	SR-BI Antibody
Catalog Number:	C48287
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SR-BI
Concentration:	1 mg/ml
Gene Name:	SCARB1 CD36L1 CLA1
Protein Name:	Scavenger receptor class B member 1 (SRB1) (CD36 and LIMPII analogous 1) (CLA-1) (CD36 antigen-like 1) (Collagen type I receptor, thrombospondin receptor-like 1) (SR-BI) (CD antigen CD36)
Other Name:	
Human Swiss Prot No.:	Q8WTV0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SR-BI at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year